

Hi Flow Stroller[®]

Liquid Oxygen Portable

A Safer Way to Breathe

Developed for oxygen patient transport when ambulation is needed, CAIRE, Inc.'s Hi Flow Stroller[®] is truly state-of-the-art cryogenics. This intrahospital system was designed to complement the existing hospital ambulatory oxygen source. The Hi Flow Stroller is considerably easier to use and, in turn, recognized as safer for both the patient and hospital staff. As a practical, money-saving solution, the effectiveness and simplicity of the Hi Flow is the answer to portable oxygen needs.

A perfect choice for hospitals and patients requiring higher flow rates, the Hi Flow Stroller is a natural fit for those patients who need a portable oxygen unit that can keep up with them. The higher oxygen flow capacity makes the Hi Flow Stroller a natural fit for activities like hiking and fishing.

Features

- Ultra-lightweight: 8 lbs. (3.6 kg)
- 1.25 L capacity
- 20 psi operating pressure
- Available in comfortable black or grey soft case
- LED meter easy to read
- Side Fill and Top Fill configurations
- Increased ease of mobility



Call 1.800.482.2473 to order!
www.CAIREmedical.com

Hi Flow Stroller®

Liquid Oxygen Portable



Product Specifications

	Hi Flow Stroller
LOX capacity	3.0 lb (1,36 kg)
Gaseous equivalent capacity	1,025 L
Weight, filled	8.5 lb (3,86 kg)
Weight, empty	5.5 lb (2,49 kg)
Height	13.5" (343 mm)
Typical use time at 2 LPM	8 hours
Operating pressure	20 psi (137 kPa)
Normal evaporation rate	1.3 lb/day (0,57 kg/day)
Standard flow control range	Continuous Flow: Off, .5, 1, 2, 2.5, 3, 4, 6, 8,10,12, 15

Ordering Information

Description	Catalog Number
Hi Flow SF (side fill) grey case	14190651
Hi Flow TF (top fill) grey case	14081227
Hi Flow SF (side fill) black case	14779019
Hi Flow TF (top fill) black case	14767421
Accessories	
Hi Flow Stroller Cart	20773616

CAIRE, Inc.
2200 Airport Industrial Drive, Suite 500
Ball Ground, GA 30107
Ph 770-721-7759 • Toll-Free 1-800-482-2473
Fax 888-932-2473

Chart BioMedical, Ltd.
Unit 6, Ashville Way
Wokingham RG41 2PL, United Kingdom
Ph +44(0) 1189 367060 • Fax +44(0) 1189 799245